# Carbon Fiber Tube (Pylon)

## Ø25/Ø30/Ø35 mm

### Instructions for Use

330126—Ø25 x 254 mm    330120—Ø30 x 254 mm    330108—Ø35 x 254 mm  
330129—Ø25 x 180 mm    330132—Ø30 x 508 mm  
330130—Ø25 x 330 mm

EN	Instructions for Use	2
FR	Instructions d'utilisation	13
DE	Gebrauchsanweisung	24
IT	Istruzioni per l'uso	35
ES	Instrucciones de uso	46
NL	Gebruiksaanwijzing	57
PL	Instrukcje użytkowania	68
PT	Instruções de utilização	79
CS	Návod k použití	90

# Contents

Contents .....	2
1 Description and Intended Purpose .....	3
2 Safety Information .....	4
3 Construction .....	5
4 Function.....	5
5 Maintenance .....	6
6 Limitations on Use .....	7
7 Fitting Advice.....	8
8 Technical Data .....	10
9 Ordering Information .....	11

---

# 1 Description and Intended Purpose

These Instructions for Use are intended for the practitioner and user unless otherwise stated.

The term *device* is used throughout to refer to 25 mm, 30 mm and 35 mm diameter Carbon Fiber Tubes (pylons).

Please read and ensure you understand these instructions, in particular all safety information and maintenance instructions.

## Application

This device is to be used exclusively for lower limb prosthetic fittings, intended for a single user.

The tube (pylon) is used as a connecting element between two Blatchford tube clamp adapters allowing rotational and length adjustment, either in a transfemoral or transtibial system.

The tubes may be cut further for customized limb alignment and length adjustment to suit the amputee.

## Features

- Carbon fiber construction
- Choice of three lengths for the 25 mm, two lengths for the 30 mm and one length for the 35 mm diameter tube.

## Activity Level

This device is suitable for Activity Levels 1-4 (weight limits apply, see *Technical Data*). Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

### Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

### Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

### Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

### Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

## Contraindications

There are no known contraindications if used in accordance with these instructions.

## Clinical Benefit

- Provides a lightweight connection between prosthetic components with additional durability.

## 2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the limb e.g. excessive play or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

Avoid exposure to extreme heat and/or cold. Do not place near any heat source. Do not leave in direct sunshine or inside a car in hot weather.

The device is designed for prolonged submersion and suitable for immersion in fresh water only. Ensure any use of the device in water complies with the conditions given in *Limitations on Use*.

The user must not adjust or tamper with the setup of the device.

The user is advised to contact their practitioner if their condition changes.

**WARNING:** Carbon fiber dust is non-toxic but is an irritant to the skin, eyes and lungs. Always read and comply with local and national Health & Safety regulations regarding the cutting of carbon fibre.

When cutting carbon fibre:

1. Ensure working area is well ventilated and/or under extraction.
2. Wear eye protection/face shield, dust mask and protective gloves.

**WARNING:** Acetone is highly flammable and toxic. Always read and comply with the solvent safety label instructions, and national and local Health & Safety regulations.

When using acetone or other solvent cleaners:

1. Ensure working area is well ventilated and/or under extraction as necessary.
2. Wear eye protection/face shield, protective gloves. Appropriate RPE (Respiratory Protective Equipment) may be required if insufficient ventilation/extraction.


---

## 3 Construction

### Principal Part

- Tube                    25 mm diameter carbon fiber tube  
                              30 mm diameter carbon fiber tube  
                              35 mm diameter carbon fiber tube

### Component Identification

---

## 4 Function

This device provides a connecting element between two Blatchford tube clamp adapters to allow for rotational and length adjustment, either of a transfemoral or transtibial system.

The tube is intended to be cut to length by a suitably trained person to suit the requirements of the amputee. This enables a comfortable fit and correct alignment to ensure correct function when using the prosthesis.

---

## 5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider e.g. unusual noises, excessive play, loss of alignment or significant wear.

Inform the practitioner/service provider of any changes in body weight and/or activity level.

### **Cleaning**

Use a damp cloth and mild soap to clean outside surfaces. DO NOT use aggressive cleansers.

***The remaining instructions in this section are for practitioner use only.***

This maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

The following routine maintenance is to be carried out at least annually:

- Check for tube tightness. If loose, remove, clean and re-fit.
- Check for defects that could affect proper function.

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check of the device is recommended and signs of wear that may affect function should be reported to their service provider (e.g. significant wear or excessive discoloration from long term exposure to UV).

Advise the user to inform the practitioner/service provider of any changes in body weight and/or activity level.

If this device is used for extreme activity, the maintenance level and interval should be reviewed and if required advice and technical support sought to plan a new maintenance schedule dependent upon the frequency and nature of the activity. This should be determined by a local risk assessment carried out by a suitably qualified individual.

---

## 6 Limitations on Use

### Intended Life

A local risk assessment should be carried out based upon activity and usage.

### Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

### Environment

This device is waterproof to a maximum depth of 1 meter. Thoroughly rinse with fresh water after use in abrasive environments such as those that may contain sand or grit, for example, to prevent wear or damage to moving parts. Thoroughly rinse with fresh water after use in salt or chlorinated water

Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).

Suitable for submersion

## 7 Fitting Advice

*The instructions in this section are for practitioner use only.*



### Cutting to Length

The device may be cut to the required length (see *Fitting Length* in Section 8). Ensure the device is cut square and the end has sharp edges removed.

**⚠ WARNING: Carbon fiber dust is non-toxic but is an irritant to the skin, eyes and lungs. Always read and comply with local and national Health & Safety regulations regarding the cutting of carbon fiber.**

**When cutting carbon fiber:**

1. Ensure working area is well ventilated and/or under extraction.
2. Wear eye protection/face shield, dust mask and protective gloves.

### Tube Preparation

Follow clamping instructions supplied with the mating component.

In instances where the patient's gait or activity may induce excessive torsional loads, it is permissible to apply a small amount of Loctite 222 to the clamping area of the tube.



**⚠ WARNING: Acetone is highly flammable and toxic. Always read and comply with the solvent safety label instructions, and national and local Health & Safety regulations.**

**When using acetone or other solvent cleaners**

1. Ensure working area is well ventilated and/or under extraction as necessary.
2. Wear eye protection/face shield, protective gloves. Appropriate RPE (Respiratory Protective Equipment) may be required if insufficient ventilation/extraction.

## 7 Fitting Advice *(continued)*

Symptom	Solution
A recurring noise occurs between the socket and the knee, at any clamping point, at the ankle or below the socket.	The device must be secured with no movement. Check that the end of the tube has been cut squarely.
The device moves out of position.	User must not use the device until adjusted, repaired or replaced. Rotation can occur if the tube has not been lightly abraded, and the tube and inside of the clamp degreased.

## 8 Technical Data

Operating and Storage Temperature Range: -15 °C to 50 °C  
(5 °F to 122 °F)

Component Weight:

Ø25 mm x 180 mm	56 g (2 oz)
Ø25 mm x 254 mm	81 g (2.9 oz)
Ø25 mm x 330 mm	109 g (3.8 oz)
Ø30 mm x 254 mm	100 g (3.5 oz)
Ø30 mm x 508 mm	206 g (7.3 oz)
Ø35 mm x 254 mm	90 g (3.2 oz)

Activity Level: 1–4

Maximum User Weight:

Ø25 mm	60 kg (132 lb)/1–4
Ø30 mm	100 kg (220 lb)/1–4
Ø35 mm	125 kg (275 lb)/1–4



Attachment Type: Proximal/Distal—Ø25, Ø30 or Ø35 mm tube

Range of Adjustment: 360° axial rotation

Length:

Ø25 mm	180 mm/254 mm/330 mm
Ø30 mm	254 mm/508 mm
Ø35 mm	254 mm

### Fitting Length

### Storage and Handling

When storing for prolonged periods, ensure the product is free from moisture and stored at room temperature.

---

## 9 Ordering Information

Device	Part Number
Carbon Fiber Tube - $\varnothing$ 25 mm x 254 mm	330126
Carbon Fiber Tube - $\varnothing$ 25 mm x 180 mm	330129
Carbon Fiber Tube - $\varnothing$ 25 mm x 330 mm	330130
Carbon Fiber Tube - $\varnothing$ 30 mm x 254 mm	330120
Carbon Fiber Tube - $\varnothing$ 30 mm x 508 mm	330132
Carbon Fiber Tube - $\varnothing$ 35 mm x 254 mm	330108

Cutting Jig	Part Number
Carbon Fiber Tube Cutting Jig - $\varnothing$ 25 mm	941238
Carbon Fiber Tube Cutting Jig - $\varnothing$ 30 mm	941237
Carbon Fiber Tube Cutting Jig - $\varnothing$ 35 mm	941257

## Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

## CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: [www.blatchford.co.uk](http://www.blatchford.co.uk)

Medical Device

Single Patient – multiple use

## Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

## Warranty

This device is warranted for 24 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement.

## Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

## Environmental Aspects

This product is fabricated from carbon fiber and cured resin that cannot be easily recycled: please dispose of it responsibly as general waste, as per local waste handling regulations.

## Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

## Trademark Acknowledgements

Blatchford is a registered trademark of Blatchford Products Limited.

## Manufacturer's Registered Address

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**Blatchford Europe GmbH**
**Endolite India Ltd.**