Silcare Walk Cushion Liner

Silcare Walk Locking Liner

User Guide

SWTTCP22–SWTTCP40
SWTTLP22–SWTTLP40

EN	User Guide	2
SL	Uporabniški priročnik	9
BG	Наръчник на потребителя	16
HR	Korisnički priručnik	23
SK	Používateľská príručka	30
HU	Felhasználói útmutató	37
EL	Οδηγός χρήσης	44
LV	Lietotāja rokasgrāmata	51
LT	Naudotojo vadovas	58
ET	Kasutaja käsiraamat	65

Contents	2
1 Description and Intended Purpose	3
2 Safety Information	3
3 Construction	4
4 Maintenance	4
4.1 Cleaning the Device	5
4.2 Cleaning the Residual Limb.....	5
5 Limitations on Use	5
6 Donning the Device	6
7 Fitting Advice.....	6
8 Technical Data	7

1 Description and Intended Purpose

These instructions are for the user. Please keep these instructions.

The term *device* is used throughout these instructions for use to refer to Silcare Walk Cushion and Locking Liner unless otherwise stated.

Make sure that you understand all instructions for use, drawing particular attention to all maintenance and safety information sections.

Application

This device is an interface component for use only as part of a lower limb prosthesis.

Intended for single user.

The device provides a cushioned socket interface and is manufactured from biocompatible materials. It is intended for a low to moderate activity user.

Note that some high activity users with poor residual limb tissue coverage and/or high sensitivity may prefer the comfort of the softer silicone Silcare Walk liner to the firmer Silcare Active liner. Conversely, some low activity users with good residual limb tissue coverage and/or low sensitivity may prefer the security of the firmer Silcare Active liner to the softer Silcare Walk liner.

Please read these instructions carefully before fitting the device.


2 Safety Information

-  This warning symbol highlights important safety information.
-  Inspect the device for damage/deterioration before and after each use. Report anything that might affect functionality.
-  Any deterioration in residuum condition or any change in sensation should be reported to the practitioner.
-  Make sure that any damaged skin or open wounds are properly and suitably dressed to prevent direct contact with the device.
-  Wearers with sensitive skin, diabetics and vascular cases should be extra vigilant and may need to apply lubricant to sensitive areas. We recommend a routine visual check and if required the user should consult their medical practitioner.
-  For other medical conditions the user should follow the advice and recommendation of a physician or a medical practitioner regarding skin care.
-  **Do not use alcohol sprays, household cleaners, or abrasives.** These cleaning materials could damage the device and irritate skin.
-  **Do not pull or stretch the device.** Fingernails, sharp jewelry, and the locking pin can tear the device. If the device is torn, stop using it and contact a Blatchford sales representative.
-  Sockets with sharp proximal edges can damage the device.
-  Take care when handling the device to avoid possible contamination with materials such as fibreglass which will stick to the device and cause skin irritation.
-  When donning a sock, clothing and the prosthetic limb, be aware that the device can build up a static charge.
-  To avoid danger of suffocation, keep the device away from babies and children.
-  Keep the device away from direct heat sources.
-  Do not overtighten the locking pin.

3 Construction

Principal Parts

- Fabric (polyamide and lycra)
- Main Body (silicone)
- Distal Cap (silicone)
- Umbrella (nylon)

4 Maintenance

Advise users to report the following to their practitioner:

- Tears in the fabric or the silicone
- Tissue damage on the residual limb
- Changes in either body weight or activity level
- Deterioration/changes to the residual limb
- Changes in the performance of the device

4.1 Cleaning the Device

Wash the inside of the device daily to avoid any build up of bacteria.

 **Take care handling the device when it is inside out to avoid picking up dust, grit and other contamination which could cause irritation to the skin.**

 **Dry thoroughly before use.**

 **Do not tumble dry.**

Washing by Hand

1. Invert the device so that its silicone side points externally.
2. Clean the silicone and the distal cap with a solution of water and unperfumed, pH balanced soap.
3. Invert the device so that its silicone side points internally.
4. Rinse the device with clean water to remove all residues.
5. Either pat dry the device with a lint-free cloth, or leave it to air dry. Take care when drying and handling the device.

Note... Always dry the device with the silicone side of the device pointing internally. Otherwise, the device can stretch and become distorted.

Machine Wash

Suitable for machine washing at 30 °C.

 **Do not invert the device for machine washing.**

4.2 Cleaning the Residual Limb

Inspect residual limb before and after prosthetic limb use or at least daily.

 **Any deterioration in residuum condition should be reported to the practitioner.**

1. Clean the skin daily with an unperfumed, pH balanced soap.
2. Rinse the skin with clean water to remove all residues.
3. Apply lotion to dry skin, as recommended by a practitioner.

 **Make sure that any damaged skin or open wound is properly and suitably dressed to prevent direct contact with the device.**

5 Limitations on Use

Environment

Avoid exposing the device to corrosive elements such as acids, industrial detergents, bleach or chlorine. The use of creams or lotions with this device should be used with caution as these can cause the device to soften and stretch or distort.

Keep away from sharp objects (such as jewelry, fingernails).

Exclusively for use between
-15°C and 50°C (5°F to 122°F).


Suitable for shower use

6 Donning the Device

 Care must be taken during donning/doffing not to damage the device with fingernails, sharp jewelry or the locking pin.

 Do NOT pull or stretch the device.

1. Invert the device so that its silicone side points externally.
2. Align the locking pin with the long axis of the residual limb. (Locking liner only.)
3. Roll the device on the residual limb while releasing any trapped air.

If the user feels numbness, tingling or any unusual sensation while using the device remove the it until normal sensation returns.

If this persists, contact your practitioner.

7 Fitting Advice



Cushion

The device must be worn in conjunction with an air-tight sleeve.


Locking

When using a locking liner the integrated matrix should finish just below fibula head/patella tendon level.

The locking pin must always follow the long axis of the residual limb. (See diagrams.)

The pin should be secured by applying Loctite 243 to the threads of the locking pin. Torque the locking pin to 3 Nm, or tighten it with your fingers by 1/8 to 1/4 of a turn.


8 Technical Data

Principal Materials	polyamide, lycra, silicone, nylon
Shore Hardness	30-35 shore 00
Component Weight (<i>size 28</i>)	695 g (1 lb 8 oz)
Activity Level	Low to Moderate
Size Range	22-40 cm
Length (See diagram below)	435 mm
Internal Length (See diagram below)	420 mm
Matrix Length (Locking Type Only)	Approx. 10 cm
Operating and Storage Temperature Range	-15 °C to 50 °C (5 °F to 122 °F)
Distal End Attachment (Locking Type Only)	M10 Requires Locking Pin With a Shoulder*
Locking Pin Shoulder Diameter	13.5-19 mm

*Locking Pin Not Supplied

Dimensions

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk

Medical Device

Single Patient – multiple use

Warranty

The device is warranted for 6 months from original date of purchase unless otherwise stated. See Blatchford catalogue for details.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

This product is fabricated from silicone rubber and fabric that cannot be easily recycled: please dispose of it responsibly as general waste, as per local waste handling regulations.

Trademark Acknowledgments

Silcare and Blatchford are registered trademarks of Blatchford Products Limited.

Manufacturer's Registered Address

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Blatchford Europe GmbH
Endolite India Ltd.