Child's 4-Bar Knee + Ø25 mm Shin Tube

Instructions for Use

019245

EN	Instructions for Use	2
SL	Navodila za uporabo	8
BG	Инструкции за употреба	14
HR	Upute za uporabu	20
SK	Návod na použitie	26
HU	Használati útmutató	32
EL	Οδηγίες χρήσης	38
LV	Lietošanas pamācība	44
LT	Naudojimo instrukcija	50
ET	Kasutusjuhend	56

Contents

Contents	2
1 Description and Intended Purpose	3
2 Safety Information	3
3 Construction	4
4 Limitations on Use	4
5 Fitting the Shin Tube	5
6 Technical Data	6
7 Ordering Information	6

1 Description and Intended Purpose

These Instructions for Use are intended for use by the practitioner and user unless otherwise stated.

The term device is used throughout these instructions for use to refer to a Child's 4-Bar Knee with a 25 mm Shin Tube.

Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.

For intended purpose, clinical benefits, contraindications, and information on how to set up and maintain the knee unit (239133), please consult the associated Instructions for Use 938017.

938017

Application

This device is to be used exclusively as part of a lower limb prosthesis.

This device is for a single user.

The carbon fiber shin tube component serves as a strong and light weight interface between the knee and ankle components of the prosthesis.

2 Safety Information


 **This warning symbol highlights important safety information.**

-  **Beware of finger trap hazard all times. The use of cosmetics is recommended to minimise the risk of finger trap.**
-  **Any changes in the performance or function of the limb e.g. instability, double-action, restricted movement, non-smooth motion or unusual noises should be immediately reported to your service provider.**
-  **Always use a hand rail when descending stairs and at any other time if available.**
-  **Any excessive changes in the heel height after finalization of alignment may adversely affect limb function.**
-  **The device is not intended for use when immersed in water or as a shower prosthesis. If the limb comes into contact with water wipe dry immediately.**
-  **Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.**
-  **Care should be taken when carrying heavy loads as this may adversely affect the stability of the device.**
-  **The user should be advised to contact their practitioner if their condition changes.**
-  **Avoid exposure to extreme heat and/or cold.**
-  **Assembly, maintenance and repair of the device must only be carried out by a suitably qualified practitioner.**

3 Construction

Principal Parts:

- Shin Tube (carbon fiber)

4 Limitations on Use

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

Environment

Avoid exposing the product to corrosive elements such as water, acids and other liquids. Also avoid abrasive environments such as those containing sand for example as these may promote premature wear.

Exclusively for use between -10 °C and 50 °C
(14 °F to 122 °F).

Suitable for outdoor use

5 Fitting the Shin Tube





The instructions in this section are for practitioner use only.

⚠ WARNING: Carbon fiber dust is non-toxic but is an irritant to the skin, eyes and lungs. Always read and comply with local and national Health & Safety regulations regarding the cutting of carbon fiber. When cutting carbon fiber:

1. Ensure working area is well ventilated and/or under extraction.
2. Wear eye protection/face shield, dust mask and protective gloves.

Cutting the Shin Tube to Length


The shin tube should be cut to the correct length to suit the user.

6 Technical Data

Operating and Storage Temperature Range:	-15 °C to 50 °C (5 °F to 122 °F)
Component Weight (Child's Knee and Shin Tube):	450 g (16 oz)
Activity Level:	1-4
Maximum User Weight:	60 kg (130 lb)
Distal Attachment:	Ø 25 mm
	Full 360° rotation
Proximal Attachment:	Blatchford Single Bolt with 7 mm A-P shift
Build Height:	135 - 407 mm

Fitting Length

7 Ordering Information

Order Example

Item	Part Number
Child's 4-Bar Knee	019245

Parts and Tools

Item	Part Number
Shin Tube 25 mm dia x 330 mm long	330130
25 mm dia Shin Tube Cutting Jig	941238

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a practitioner.

Warranty

This device is warranted for 24 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

Where possible the components should be recycled in accordance with local waste handling regulations.

Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

Manufacturer's Registered Address

Blatchford Products Limited, Lister Road, Basingstoke RG22 4AH, UK.

Trademark Acknowledgements

Blatchford is a registered trademark of Blatchford Products Limited.

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Blatchford Europe GmbH
Endolite India Ltd.