# Silcare Breathe Active Cushion Liner

## User Guide

SBATTCP22-SBATTCP40

EN	Instructions for Use	2
FR	Instructions d'utilisation	10
DE	Gebrauchsanweisung	18
IT	Istruzioni per l'uso	26
ES	Instrucciones de uso	34
NL	Gebruiksaanwijzing	42
PL	Instrukcje użytkownika	50
PT	Instruções de utilização	58
CS	Návod k použití	66

# Contents

Contents .....	2
1 Description and Intended Purpose .....	3
2 Safety Information .....	4
3 Construction .....	5
4 Maintenance .....	5
4.1 Cleaning the Device .....	6
4.2 Cleaning the Residual Limb.....	6
5 Limitations on Use .....	6
6 Donning the Device .....	7
7 Fitting Advice.....	7
7.1 Moisture Collects Inside the Distal end of the Device .....	7
7.2 Any Discoloration is Noticed on the Residuum.....	7
7.3 Oedematous Swellings Corresponding to Distal Perforations .....	7
8 Technical Data .....	8

---

# 1 Description and Intended Purpose

These instructions are for the user. Please keep these instructions.

The term *device* is used throughout these instructions for use to refer to Silcare Breathe Active Cushion Liner.

Make sure that you understand all instructions for use, drawing particular attention to all maintenance and safety information sections.

## **Application**

This device is an interface component for use only as part of a lower limb prosthesis.

Intended for single user.

It is intended for moderate to high activity users with well covered residual limbs.

The device provides a perforated cushioned socket interface which allows moisture to pass through the perforations keeping the skin dry. It is manufactured from biocompatible materials.

Note that some high activity users with poor residual limb tissue coverage and/or high sensitivity may prefer the comfort of the softer silicone Silcare Walk liner to the firmer Silcare Active liner. Conversely, some low activity users with good residual limb tissue coverage and/or low sensitivity may prefer the security of the firmer Silcare Active liner to the softer Silcare Walk liner.

## 2 Safety Information


-  This warning symbol highlights important safety information.
-  Inspect the device for damage/deterioration before and after each use. Report anything that might affect functionality.
-  Any deterioration in residuum condition should be reported to the practitioner. If the residual limb has tissue damage, the use of the device should be stopped.
-  Make sure that any damaged skin or open wound is properly and suitably dressed to prevent direct contact with the device.
-  Wearers with sensitive skin, diabetics and vascular cases should be extra vigilant and may need to apply lubricant to sensitive areas. We recommend a routine visual check and if required the user should consult their medical practitioner.
-  For other medical conditions the user should follow the advice and recommendation of a physician or a medical practitioner regarding skin care.
-  Enlarged perforations can trap the skin and cause blisters. If the perforations enlarge, stop using the device.
-  If distal oedematous swellings corresponding to distal perforations in the liner occur, use of the liner should be discontinued and the swellings reported to the practitioner.
-  **Do not** use alcohol sprays, household cleaners, or abrasives. These cleaning materials could damage the device and irritate skin.
-  **Do not** pull or stretch the fabric. Fingernails and sharp jewelry, can tear the fabric. If the fabric is torn, stop using it and contact a Blatchford sales representative.
-  Sockets with sharp proximal edges can tear the device.
-  Take care when handling the device to avoid possible contamination with materials such as fibreglass which will stick to the device and cause skin irritation.
-  When donning a sock, clothing and the prosthetic limb, be aware that the device can build up a static charge.
-  To avoid danger of suffocation, keep the device away from babies and children.
-  Keep the device away from direct heat sources.
-  **Only** use the device in combination with corrosion resistant components.

---

## 3 Construction

### Principal Parts

- Fabric (polyamide and lycra)
- Main Body (silicone)
- Distal Cap (silicone)


---

## 4 Maintenance

Advise users to report the following to their practitioner:

- Tears in the fabric or the silicone
- Tissue damage on the residual limb
- Distal oedematous swellings
- Enlarged perforations
- Changes in either body weight or activity level
- Deterioration/changes to the residual limb
- Changes in the performance of the device
  - Moisture in the distal end of the device

Note... Perspiration discolors some socket materials.

## 4.1 Cleaning the Device

Wash the inside of the device daily to avoid any build up of bacteria.

 **Take care handling the device when it is inside out to avoid picking up dust, grit and other contamination which could cause irritation to the skin.**

 **Dry thoroughly before use.**

 **Do not tumble dry.**

### Washing by Hand

1. Invert the device so that its silicone side points externally.
2. Clean the silicone and the distal cap with a solution of water and unperfumed, pH balanced soap.
3. Invert the device so that its silicone side points internally.
4. Fill the device with warm water while both holding the proximal end closed and squeezing the distal end to flush the water through the perforations.
5. Rinse the device with clean water to remove all residues.
6. Either pat dry the device with a lint-free cloth, or leave it to air dry. Take care when drying and handling the device.

Note... Always dry the device with the silicone side of the device pointing internally. Otherwise, the fabric can stretch and become distorted.

### Machine Wash

Suitable for machine washing at 30 °C.

 **Do not invert the device for machine washing.**

## 4.2 Cleaning the Residual Limb

Inspect residual limb before and after prosthetic limb use or at least daily.

 **Any deterioration in residuum condition should be reported to the practitioner. If the residual limb has tissue damage, the use of the device should be stopped.**

1. Clean the skin daily with an unperfumed, pH balanced soap.
2. Rinse the skin with clean water to remove all residues.
3. Dry residual limb.
4. Apply lotion to dry skin, as recommended by a practitioner.

 **Make sure that any damaged skin or open wound is properly and suitably dressed to prevent direct contact with the device.**

---

## 5 Limitations on Use

### Intended Life

A local risk assessment should be carried out based upon activity and usage.

### Environment

Avoid exposing the device to corrosive elements such as acids, industrial detergents, bleach or chlorine. The use of creams or lotions with this device should be used with caution as these can cause the device to soften and overstretch or distort.

Keep away from sharp objects (such as jewelry, fingernails).

Exclusively for use between -15 °C and 50 °C (5 °F to 122 °F).

Suitable for shower use

---


## 6 Donning the Device

 **Care must be taken during donning/doffing not to damage the device with fingernails or sharp jewelry.**

 **Do NOT pull or stretch the device.**

1. Invert the device so that its silicone side points externally.
2. Roll the device on the residual limb while releasing any trapped air.

After donning the device, allow the user to wear it for 10 minutes. If the user feels numbness, tingling or any unusual sensation within 10 minutes, do the following:

Note... Doff the device, and wait until normal sensation returns.

3. Don the device again.
4. If the user feels numbness, tingling or any unusual sensation again, stop using the device.
5. Add air-tight sleeve of choice.

---

## 7 Fitting Advice

### 7.1 Moisture Collects Inside the Distal end of the Device

Cause	Solution
The perforations are clogged.	Clean the device. (See <i>Cleaning the Device</i> section)

### 7.2 Any Discoloration is Noticed on the Residuum

Cause	Solution
Incorrectly fitting device.	Contact your practitioner.


### 7.3 Oedematous Swellings Corresponding to Distal Perforations

Cause	Solution
Excessive distal contact.	Reduce end pressure by addition of extra socks or lengthen/remake socket.

## 8 Technical Data

Principal Materials	polyamide and lycra, silicone
Shore Hardness	40 Shore 00
Component Weight (size 28)	695 g (1 lb 8 oz)
Activity Level	3-4
Size Range	22-40 cm
Length (See diagram below)	435 mm
Internal Length (See diagram below)	420 mm
Operating and Storage Temperature Range	-15 °C to 50 °C (5 °F to 122 °F)

### Dimensions

## Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

## CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: **[www.blatchford.co.uk](http://www.blatchford.co.uk)**.

Medical Device

Single Patient – multiple use

## Warranty

The device is warranted for 6 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses, and exemptions.

See Blatchford website for the current full warranty statement.

## Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

## Environmental Aspects

This product is fabricated from silicone rubber and fabric that cannot be easily recycled: please dispose of it responsibly as general waste, as per local waste handling regulations.

## Trademark Acknowledgments

Silcare Breathe and Blatchford are registered trademarks of Blatchford Products Limited.

## Manufacturer's Registered Address

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© Blatchford Products Limited 2023. All rights reserved.
**Endolite India Ltd.**
Blatchford Europe GmbH