# KX06

## Instructions for Use

KX06V2

|    |                       |     |
|----|-----------------------|-----|
| EN | Instructions for Use  | 2   |
| DA | Brugsanvisning        | 19  |
| NO | Bruksanvisning        | 36  |
| FI | Käyttöohjeet          | 53  |
| SV | Bruksanvisning        | 70  |
| RU | Инструкция протезиста | 87  |
| ZH | 使用说明                  | 104 |

|  |    |
|--|----|
| Contents .....                                     | 2  |
| Package Contents .....                             | 2  |
| 1 Description and Intended Purpose .....           | 3  |
| 2 Safety Information .....                         | 4  |
| 3 Construction .....                               | 5  |
| 4 Function.....                                    | 6  |
| 5 Maintenance.....                                 | 6  |
| 6 Limitations on Use .....                         | 7  |
| 7 Bench Alignment.....                             | 8  |
| 8 Static Alignment .....                           | 9  |
| 8.1 Check Anterior-Posterior Alignment.....        | 9  |
| 8.2 Adjusting Yield/Stance Resistance.....         | 9  |
| 9 Dynamic Alignment: Swing Resistance.....         | 10 |
| 9.1 Flexion Resistance.....                        | 10 |
| 9.2 Extension Resistance .....                     | 11 |
| 9.3 Stirrup Function and Features .....            | 11 |
| 9.4 Static Alignment Check .....                   | 12 |
| 9.5 Dynamic Alignment Adjustment Flow Diagram..... | 13 |
| 10 Fitting Advice.....                             | 14 |
| 11 Socket Attachment.....                          | 14 |
| 12 Fitting the Cosmetic Cover .....                | 15 |
| 13 Knee Pad Removal/Replacement .....              | 15 |
| 14 Fitting T-Slot Style Adapters .....             | 16 |
| 15 Technical Data .....                            | 16 |
| 16 Ordering Information .....                      | 17 |

## Package Contents

1. KX06V2 Knee prosthesis
2. Cosmetic cover
3. Clinician's manual
4. User guide
5. Parts bag with:
  - i) 5 x M4 x 12 mm low profile hex socket cap head screws
  - ii) 5 x Cover fixing spacers
  - iii) Hex key, 2.5 mm

---

# 1 Description and Intended Purpose

These instructions are for the practitioner.

The term *device* is used throughout these instructions for use to refer to KX06V2.

## Application

This device is a knee unit that is for use exclusively as a component of a lower limb prosthesis.

This device is a polycentric knee unit that uses both a 4-bar mechanism and a hydraulic cylinder to control the phases of the gait cycle. It provides swing and stance control for active hip disarticulation, transfemoral and knee disarticulation users.

The geometry of the device provides increased toe clearance during swing phase. The build height is short and upon sitting minimizes protrusion of the prosthetic knee beyond the contralateral limb. This combination of features reduces asymmetry for users with longer residual limbs.

This device is manufactured from water resistant materials and is designed to accommodate limited exposure to such environments.

## Features

- Adjustable stance support
- Adjustable swing (flexion and extension) resistances
- Reduced build height\*
- Increased toe clearance during swing phase\*
- Reduced protrusion during sitting\*
- Flexion Lock mode
- Cycling mode

\*Compared to a single axis knee unit.

## Clinical Benefits

- Improved ground clearance reducing trips and falls risk
- Geometry increases knee stability during stance phase

## Activity Level

This device is for Activity Level 3 and 4 users.

## Contraindications

This device might not be suitable for Activity Level 1–2 users or for use in competitive sports events. These types of users might be better served by a specially designed prosthesis that is optimized for their needs.

This device is not suitable for users with:

- Residual muscular weakness, contractures that cannot be correctly accommodated, or proprioceptive dysfunction including poor balance
- Contralateral joint instabilities or pathology
- Complicated conditions involving multiple disabilities

This device is for a single user.

Users must be given gait training before using this device.

Make sure that the user has understood all instructions for use, drawing particular attention to all maintenance and safety information sections.

---

## 2 Safety Information

**This warning symbol highlights important safety information which must be followed carefully.**

**1. Be aware of finger trap hazard at all times.**

**2. Any changes in the performance or function of the limb e.g. instability, double-action, restricted movement, non-smooth motion or unusual noises should be immediately reported to your service provider.**

**3. Always use a hand rail when descending stairs and at any other time if available.**

**4. Any excessive changes in heel height after finalization of alignment may adversely affect limb function.**

**5. After continuous use the cylinder can become hot.**

**6. The device is not designed for prolonged submersion. Ensure any use of the device in water complies with the conditions given in "Limitations on Use" (see Section 6).**

**7. After any period of immersion in water, rinse the device with fresh clean water and then thoroughly dry before reuse.**

**8. Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.**

**9. Care should be taken when carrying heavy loads as this may adversely affect the stability of the device.**

**10. The user should be advised to contact their practitioner if their condition changes.**

**11. Avoid exposure to extreme heat and/or cold.**

**12. Assembly, maintenance and repair of the device must only be carried out by a suitably qualified clinician.**


Note... The device is supplied with the cosmetic cover detached. Although the device can be configured with the cosmetic cover attached, the adjustment cap is more accessible with the cover detached. Therefore it is easier to first align and configure the device for the user (see Sections 7 to 10) while the cosmetic cover is off (taking care to avoid finger trap hazards), then, when configuration is completed, immediately attach the cosmetic cover (see Section 12) to help protect against finger trap hazards in use.

# 3 Construction


## Principal Parts

- Frame Aluminium alloy, brass, stainless steel, titanium, steel
- Knee control Various materials, principally aluminium alloy, stainless steel, acetal homopolymer, nitrile butadiene rubber (NBR), hydraulic fluid
- Knee pad Polyurethane
- Cosmetic cover Thermoplastic elastomer

## Component Identification

## Rear View Knee Control Detail

## Cosmesis

The following standard Blatchford continuous polyester foam cosmesis may be used with this device:

561021 - Fairing Unshaped Long

---

## 4 Function

The knee is a polycentric design comprising pivots connecting a chassis to a carrier via four link arms. A hydraulic knee control sits within this frame assembly, fitted between the chassis posterior pivot (rear proximal axis) and the carrier anterior pivot (distal axis).

The device provides increased toe clearance during swing phase as a result of ease of initiation of free swing and mechanical advantage of the geometry compared to a single axis device, and short protrusion when sitting.

The device has dual stance control, combining the geometric stability of the polycentric knee linkage with an adjustable hydraulic yielding support activated from knee extension prior to heel strike and released by a hyper-extension moment at any time, usually occurring after mid-stance.

The device provides an adjustable hydraulic swing phase control to accommodate variable cadence.

---

## 5 Maintenance

An annual visual inspection is recommended. Check for visual defects that may affect proper function.

Maintenance must be carried out by competent personnel.

Do not carry out maintenance on this device yourself. Instead, return it to a Blatchford center for maintenance. If the device is still under warranty, we will loan you another knee unit while we carry out maintenance.

The user should be advised:

Any changes in performance of this device must be reported to the practitioner.

Changes in performance may include:

- Increase or decrease in knee stiffness
- Instability
- Reduced knee support (free movement)
- Any unusual noise\*

\* Due to the nature of the hydraulics there may be some slight air noise from the cylinder during the first few steps. This is not detrimental to the function of the unit and should dissipate quickly. If symptoms persist, please consult your practitioner.

### **Storage and Handling**

When storing for prolonged periods place the device vertically with trunnion uppermost. Alternatively and in addition, air management may be assisted by cycling the knee control several times with the stance resistance switched off (refer to Section 9.3).

Use product packaging supplied.

### **Cleaning**

Use a damp cloth and mild soap to clean the outside surfaces.

DO NOT use aggressive cleansers.

---

## 6 Limitations on Use

### Intended Life

A local risk assessment should be carried out based upon activity and usage.

### Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

### Environment

The device is suitable for submersion in fresh water only.

Knee must be locked while submerged, then rinsed and thoroughly dried.

Light surface corrosion affects neither the function nor the security of this device.

However, if it is heavily corroded, stop using and contact your practitioner.

Avoid exposing the device to abrasive or corrosive environments, such as those containing sand, salt or chemically treated water e.g. swimming pools. After any exposure to such environments rinse the device in fresh water and dry thoroughly.

Exclusively for use between -10 °C and 50 °C (-14 °F and 122 °F).


Suitable for submersion

It is recommended that only Blatchford products be used in conjunction with the device.

# 7 Bench Alignment

Users be aware of potential finger trap hazard.

## Additional Components

Follow the alignment instructions that are supplied with any additional components.

# 8 Static Alignment


To achieve optimal function from the Blatchford hydraulic knee control the knee must be aligned geometrically **stable**.

Check flexion is fully accommodated when worn by user.

## 8.1 Check Anterior-Posterior Alignment


**1**

- 1) Ensure stirrup is down.
- 2) Stand with feet side by side holding onto rails.
- 3) Load prosthesis.

**Knee *stable* during weight bearing:**


- 1) Proceed to adjusting yield/stance resistance, see Section 8.2.

OR

**Knee *unstable* during weight bearing:**

- 1) Check flexion is fully accommodated.
- 2) Confirm bench alignment, see Section 7.
- 3) Slide socket anteriorly to make knee just stable.


**Stance Resistance (Yield)**

**Do not turn the yield adjustment screw counterclockwise beyond 180° position as this will cause the control unit to operate inconsistently (may or may not release from stance.)**

## 8.2 Adjusting Yield/Stance Resistance


**2**


2 mm Hex Key

Stirrup

**Yield Adjuster**

9 = Maximum Stance Resistance


1 = Minimum Stance Resistance

**2**

Having first checked A-P alignment as above (initial adjustment position is factory set):


- 1) Ensure stirrup in **down** position.
- 2) Standing with feet side by side try to sit (do not hyperextend limb).

If excessive resistance, reduce resistance by turning the yield adjustment screw with the hex key counterclockwise until correct resistance is felt to sitting naturally.

# 9 Dynamic Alignment: Swing Resistance

## 9.1 Flexion Resistance



Raise latch to allow adjustment cap to turn.

Lower latch to lock adjustment cap.

Note... If you cannot turn the adjustment cap with your fingers, use a hydraulic-knee-control adjuster (940091).


**!** Take care not to force the adjuster ring and overcome the adjustment stops; this could cause loss of function.

Note... The same adjustment cap is used to adjust both flexion resistance and extension resistance.

### 3


- 1) Initial adjustments are factory set (Extension 2, Flexion 4).
- 2) Observe the user walking.
- 3) If there is excessive heel rise: Increase flexion resistance.
- 4) If there is insufficient heel rise: Decrease flexion resistance.

**!** **Warning!** If adjusting the flexion resistance has no effect on heel rise, check stirrup is in down position and ensure that the user is initiating flexion whilst still loading the toe. Note that a hyperextension moment is needed about the knee to initiate flexion.

The flexion scale on the adjustment cap is numbered from left to right 1 (minimum) to 10 (maximum). Directly below the flexion adjustment cap is an indicator mark.

Turning the cap counterclockwise reduces flexion resistance until the minimum resistance is reached where "1" on the flexion scale is over the indicator mark.

Indicator Mark  
Minimum Flexion (1)      Maximum Flexion (10)

Turning the cap clockwise increases the resistance to a maximum setting at "10".

To adjust flexion resistance:


- 1) Raise the latch.
- 2) Turn the adjustment cap to the required flexion setting.
- 3) Lower the latch to lock the cap and prevent accidental adjustment.

## 9.2 Extension Resistance

**4**

- 1) Observe the user walking.
- 2) If there is excessive terminal impact on knee extension: **increase** resistance.
- 3) If the knee does not extend satisfactorily: **decrease** resistance.

Note... As a 'rule of thumb' flexion resistance should be greater than extension resistance.


Note... If there is any doubt over the position of the indicator ring in relation to the settings, its position can be reset by turning the adjustment cap to maximum flexion, then maximum extension prior to making any adjustments.

After completing all flexion and extension resistance adjustments, refit the cosmetic cover if previously removed (see Section 12).

**Make a note of the flexion setting before starting to adjust the extension resistance.**

To increase extension resistance:

- 1) Raise the latch and turn the adjustment cap in the clockwise direction until it reaches 10 on the flexion scale.

- 2) Continue turning the cap clockwise to engage and turn the indicator ring clockwise; this moves the indicator mark from right to left (clockwise), increasing extension resistance from 1 to 10 on the extension scale.

Indicator Mark  
Minimum Extension (1)      Maximum Extension (10)

To reduce extension resistance:

- 3) Raise the latch and turn the adjustment cap anticlockwise to minimum on the flexion scale then continue turning until the desired extension resistance is reached on the extension scale.


When the correct extension resistance setting is achieved, turn the adjustment cap back to the previously noted flexion resistance setting, taking care not to re-engage and turn the indicator ring again.

Lower the latch to lock the cap and prevent accidental adjustment.

## 9.3 Stirrup Function and Features


**How to switch off stance resistance (for cycling etc.)**

- 1) Load toe of prosthesis.
- 2) While toe is loaded, lift stirrup up.
- 3) Stance now off but knee may be unstable.
- 4) Lower stirrup to resume normal operation.


**How to lock knee against flexion (for prolonged standing etc.)**


- 1) Flex knee slightly without overriding stance resistance.
- 2) Lift stirrup.
- 3) The knee is now locked against flexion but will extend.
- 4) Lower stirrup to resume normal operation.

## 9.4 Static Alignment Check

## 9.5 Dynamic Alignment Adjustment Flow Diagram


## 10 Fitting Advice

| Problem  | Solution  |
|--|---|
| The user's sitting movement is too slow.       | Decrease the stance support (see Section 8.2).  |
| The user's sitting movement is too fast.       | Increase the stance support (see Section 8.2).  |
| The device does not release from stance.       | <ol style="list-style-type: none"><li>1. Make sure that the user loads the toe of the foot unit while initiating swing and that the device reaches full extension.</li><li>2. Make sure that the ankle-foot unit's dorsiflexion resistance is not too low.</li><li>3. Re-check alignment. Socket may be too posteriorly aligned causing the knee to flex/yield at heel strike).</li></ol> |
| The heel rise is too high during swing.        | Increase the flexion resistance (see Section 9.1).  |
| The heel rise is too low during swing.         | Decrease the flexion resistance (see Section 9.1).  |
| The user is experiencing terminal impact.      | Increase the extension resistance (see Section 9.2).  |
| The device does not fully extend during swing. | Decrease the extension resistance (see Section 9.2).  |

## 11 Socket Attachment

Ensure socket does not contact the knee control or posterior link arms when fully flexed.  
The socket may contact the carrier only.


# 12 Fitting the Cosmetic Cover

The cover and mounting screws are supplied detached allowing easy access to the alignment controls.

**!** Always fit the cosmetic cover after all alignment adjustments have been made to help prevent finger trap hazards.

1. Fit the cosmetic cover to the device at the chassis anterior pivots and carrier anterior pivots using the 4 x M4 hex screws and 4 x spacers supplied (see figure below).
2. Tighten the screws using the 2.5 mm hex key supplied. Avoid over-tightening the screws.

**!** Avoid over-tightening the cover screws

Fitting Cosmetic Cover

# 13 Knee Pad Removal/Replacement


## Knee Pad Removal

1. Using the hex key supplied, unscrew, remove and set aside the two lateral M3 hex screws that retain the knee pad.
2. Slide the knee pad horizontally off the chassis.

## Knee Pad Replacement

1. Slide the replacement knee pad horizontally onto the chassis until the lateral knee pad and chassis screw holes align (see figure opposite).
2. Using the hex key supplied, fasten the knee in place with the two M3 hex screws. Avoid over-tightening the screws.

**!** Avoid over-tightening the knee pad screws

M3 Hex Screw  
(x2)  
910408


Knee Pad Replacement

# 14 Fitting T-Slot Style Adapters

If fitting a T-slot style adapter:

1. Remove knee pad as described in Section 13.
2. Unscrew two M3 x 8 hex socket set screws (see figure right) until the T-Slot opening is unobstructed.
3. Slide adapter into T-slot.
4. Screw in the two M3 x 8 screws until tight against the chassis to prevent the adapter from sliding out of the T-slot. Do not over-tighten screws.
5. Refit knee pad as described in Section 13.



 **Avoid over-tightening the screws after fitting T-slot adapter**

**Unscrewing T-Slot Screws Before Fitting Adapter**

# 15 Technical Data

|  |  |
|--|--|
| Operating and Storage Temperature Range: | -10 °C to 50 °C<br>(14 °F to 122 °F)               |
| Component Weight:                        | 1.34 kg (2lb 15 oz)                                |
| Activity Level:                          | 3-4  |
| Maximum User Weight:                     | Level 3: 150 kg (330lb)<br>Level 4: 145 kg (320lb) |
| Maximum Flexion Angle:                   | 160 °  |
| Proximal Alignment Attachment:           | 4-bolt interface<br>T-slot interface               |
| Distal Alignment Attachment:             | Blatchford male pyramid                            |
| Build Height (see Diagram):              | 220 mm (8.7 inches)                                |

## 16 Ordering Information

| Description   | Part Number |
|---|-------------|
| <b>Device</b>   |             |
| KX06V2  | KX06V2      |
| <b>Spare Parts</b>  |             |
| Knee pad kit (includes knee pad, 2.5mm hex key and 2 x M3 x 10 mm hex socket cap head screws)   | 239561      |
| Knee pad screw (M3 x10mm hex socket cap head screw)   | 910408      |
| Cover replacement kit (includes cosmetic cover, 2.5 mm hex key, 5 x M4 x 12 mm low profile hex socket cap head screws and 5 x cover fixing spacers) | 239661      |
| Spacer and screws kit for cosmetic cover (includes 4 x M4 x 12 mm low profile hex socket cap head screws and 4 x cover fixing spacers)              | 239761      |
| Cover screw: M4 x 12 mm low profile hex socket cap head screw   | 913216      |
| Cover spacer  | 236122      |
| Hydraulic-knee-control adjuster   | 940091      |

### Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

### CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: [www.blatchford.co.uk](http://www.blatchford.co.uk)

Medical Device

Single Patient – multiple use

### Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

## **Warranty**

This device is warranted for 36 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement.

## **Reporting of Serious Incidents**

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

## **Environmental Aspects**

This product contains hydraulic oil, mixed metals and plastics. It should be recycled where possible in accordance with local waste recycling regulations.

## **Retaining the Packaging Label**

You are advised to keep the packaging label as a record of the device supplied.

## **Trademark Acknowledgements**

KX06 and Blatchford are registered trademarks of Blatchford Products Limited.

## **Manufacturer's Registered Address**

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**Blatchford Europe GmbH**
**Endolite India Ltd.**