Multiflex Standard

519119-519136

Multiflex Slim

519141-519148

Instructions for Use

EN	Instructions for Use	2
SL	Navodila za uporabo	12
BG	Инструкции за употреба	22
HR	Upute za uporabu	32
SK	Návod na použitie	42
HU	Használati útmutató	52
EL	Οδηγίες χρήσης	62
LV	Lietošanas pamācība	72
LT	Naudojimo instrukcija	82
ET	Kasutusjuhend	92

Contents

Contents	2
1 Description and Intended Purpose	3
2 Safety Information	4
3 Construction	5
4 Function.....	6
5 Maintenance	6
6 Limitations on Use	7
7 Bench Alignment.....	8
7.1 Static Alignment.....	8
7.2 Dynamic Alignment	8
8 Technical Data	9
9 Ordering Information	10

1 Description and Intended Purpose

These Instructions for Use are intended for use by the practitioner and user unless otherwise stated.

The term *device* is used throughout these instructions for use to refer to Multiflex feet.

Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.

Application

This device is to be used exclusively as part of a lower limb prosthesis.

Intended for a single user.

This combination provides the user with multi-axial motion when used in conjunction with a multiflex ankle, which enhances stability. The low profile design with sandal toe gives the foot a more natural appearance.

Activity Level

This device is recommended for Activity Level 1–2 users.

There are exceptions to our recommendations and we want to allow for individual circumstances. There may be Activity Level 3–4* users who may benefit from this product, but this decision should be made with sound and thorough justification.

*Maximum user weight: 100 kg (220 lb).

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low level environmental barriers such as, curbs, stairs or uneven surfaces. Typical of the limited community ambulator.

Activity Level 3

Has the ability or potential for ambulation with variable cadence.

Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

Contraindications

This device may not be suitable for use in competitive sports events. This type of user may be better served by a specially designed prosthesis optimized for their needs.

Clinical Benefits

- Low stiffness at weight acceptance leads to early foot-flat and greater stability for lower mobility users, compared to single axis feet and Solid-Ankle-Cushioned-Heel feet
- No loss of stability for bilateral users, compared to rigid ankles
- Easier to walk on uneven ground for bilateral users, compared to rigid ankles
- Easier to walk on slopes for bilateral users, compared to rigid ankles
- Easier to walk longer distances for bilateral users, compared to rigid ankles
- Improved stance phase timing symmetry, compared to rigid ankles
- Reduced sound limb loading, compared to rigid ankles

2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the limb e.g. restricted movement, excessive movement, non-smooth motion or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

The device is not suitable for extreme sports, running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk. Recreational cycling is acceptable.

Assembly, maintenance and repair of the device must only be carried out by a practitioner or suitably trained technician.

The user should contact their practitioner if their condition changes.

Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.


To minimise the risk of slipping and tripping, appropriate footwear that fits securely onto the footshell must be used at all times.

3 Construction

Principal Parts

- Foot Shell - Polyurethane
- Molded Keel - Fiber reinforced thermoplastic/aluminium alloy
- Foot Bolt - Stainless steel
- Concave Washer - Stainless steel

Component Identification

4 Function

The device combines with a resilient multi-axial Multiflex ankle to give a smooth, stable walking action and allow some compliance to uneven surfaces. The device comprises a lightweight, fiber reinforced thermoplastic keel which gives a gentle rolling action through stance and features a sandal toe for a more natural appearance—ideal for lower level activity users.

A range of heel heights can be accommodated via the serrated ankle interface. Multiflex Slim features additional heel height capacity for a higher heel.

5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider e.g. unusual noises, increased stiffness or restricted/excessive movement, significant wear or excessive discoloration from long term exposure to UV.

Inform the practitioner/service provider of any changes in body weight and/or activity level.

Cleaning

Use a damp cloth and mild soap to clean outside surfaces, do not use aggressive cleansers.

The remaining instructions in this section are for practitioner use only.

This maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

It is recommended that the following maintenance is carried out annually:

- Check the foot bolt for tightness

Advise users to report the following to their medical practitioner:

- Changes in body weight or activity level
- Changes in the performance of this device, for example unusual noises or play

Make sure the user has understood all instructions for use, drawing particular attention to the Maintenance section.

6 Limitations on Use

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Lifting Loads

User weight and activity is governed by the stated limits. Load carrying by the user should be based on a local risk assessment.

Environment

Avoid exposing the Multiflex foot to corrosive elements such as water, acids and other liquids. Avoid abrasive environments such as those containing sand, as these may promote premature wear.


Exclusively for use between -10 °C and 50 °C
(14 °F to 122 °F).

Suitable for outdoor use

7 Bench Alignment

The instructions in this section are for practitioner use only.

Align trans-femoral devices according to fitting instructions supplied with the knee, keeping the build line relative to the Multiflex ankle, as shown.

7.1 Static Alignment

With the flexion accommodated the build line should pass through the center of the socket in the sagittal (A-P) plane and through the Multiflex foot, as shown. Heel height of the shoe should be accommodated.

7.2 Dynamic Alignment

Coronal Plane

Make sure that M-L thrust is minimal by adjusting relative positions of socket and foot.

Sagittal Plane

Check for smooth transition from heel strike to toe-off. Also, make sure that when standing, the heel and toe are evenly loaded and that both are touching the floor.

Check that the serrations are engaged before tightening the foot bolt.

Incorrectly engaged serrations can lead to noise and result in sudden loosening of the foot. Confirm correct engagement of serrations by slackening foot bolt 1/8 turn. Foot should be loose but not rotate. then torque to correct torque.


* Allow for the user's own footwear

8 Technical Data


Operating and Storage Temperature Range:	-15 °C to 50 °C (5 °F to 122 °F)
Component Weight (size 26):	375 g (13.2 oz)
Activity Level:	1–2
Maximum User Weight:	
Activity Level 1-2	125 kg (275 lb)
Activity Level 3-4	100 kg (220 lb)
Build Height:	40 mm
Nominal/Mid-range Heel Height:	Multiflex Slim — 25 mm Multiflex Standard — 20 mm
Range of Heel Height Adjustment:	± 15 mm (3x5 mm increments)

Fitting Length

Multiflex Standard

Multiflex Slim

9 Ordering Information

Multiflex Standard

Size (cm)	Part Number — Left	Part Number — Right
22	519119	519120
23	519121	519122
24	519123	519124
25	519125	519126
26	519127	519128
27	519129	519130
28	519131	519132
29	519133	519134
30	519135	519136

Multiflex Slim

Size (cm)	Part Number — Left	Part Number — Right
22	519141	519142
23	519143	519144
24	519145	519146
25	519147	519148

Part	Part Number
Foot Bolt	400022
Concave Washer	500021

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

Warranty

This device is warranted for 24 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

The above may vary by market; consult your local representative for details.

See the Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

To help prevent potential harm to the environment or to human health from uncontrolled waste disposal, Blatchford offers a take-back service. Please contact Customer Services for details.

Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgements

Multiflex and Blatchford are registered trademarks of Blatchford Products Limited.

Manufacturer's Registered Address

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Blatchford Europe GmbH
Endolite India Ltd.