# BK Top Housing

## Instructions for Use

339008 (30 mm)

EN	Instructions for Use	2
DA	Brugsanvisning	13
NO	Bruksanvisning	24
FI	Käyttöohjeet	35
SV	Bruksanvisning	46
RU	Инструкция протезиста	57
ZH	使用说明	68
AR	إرشادات الاستخدام	79

# Contents

Contents .....	2
1 Description and Intended Purpose .....	3
2 Safety Information .....	4
3 Construction .....	5
4 Function.....	6
5 Maintenance .....	6
6 Limitations on Use .....	7
7 Bench Alignment.....	8
8 Fitting Advice.....	9
9 Technical Data .....	10
10 Ordering Information .....	11

---

# 1 Description and Intended Purpose

These Instructions for Use are intended for use by the practitioner and user unless otherwise stated.

The term *device* is used throughout these instructions for use to refer to the BK Top Housing (BK=Below Knee).

Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.

## Application

The Device is designed for use exclusively as part of a lower limb prosthesis.

This device is a 30 mm tube clamp connector allowing connection of Blatchford single bolt devices to a 30 mm tube.

Once the secondary device is attached there will be facility for:

- $\pm 5$  mm horizontal shift
- Rotation (up to 360°)

Intended for single user only.

## Activity Level

This device is suitable for activity levels 1 to 3; weight limits apply, see Technical Data, Section 9. Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

### Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited ambulator.

### Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

### Activity Level 3

Has the ability or potential for ambulation with variable cadence.

Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

## Contraindications

This device is not suitable for immersion.

Maximum user weights apply, see Technical Data Section.

## Clinical Benefits

- Allows connection of single bolt devices to 30mm tube

## 2 Safety Information

**This warning symbol highlights important safety information which must be followed carefully.**

Any changes in the performance or function of the device e.g. unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

Avoid exposure to extreme heat and/or cold.

To reduce the risk of injury due to failure or loosening of the screw connections ensure the bolt threads are cleaned thoroughly before each installation.

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified clinician that has attended an approved training course.

The user must not adjust or tamper with the setup of the device.

The user should be advised to contact their practitioner if their condition changes.

Always apply the specified torque values to the screws/ nuts. Never replace the screws with alternative screws. Always use the bolts specified.

Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.


Make sure that no liquids enter the device during donning or normal use.

# 3 Construction

## Principal Parts

- Tube Clamp                      Aluminium
- M6 Screw                         Carbon Steel
- M6 Lock washer                 Steel

## Component Identification

---

## 4 Function

The device provides a connection point between a single bolt socket adapter and a 30 mm tube.

---

## 5 Maintenance

Visually check the device regularly.

Report any changes in performance and signs of wear that may affect function of this device to the practitioner/service provider.

Changes in performance may include:

- Instability
- Any unusual noise
- Excessive play or loss alignment

Signs of wear that may affect function may include:

- Surface damage
- Corrosion

Inform the practitioner/service provider of any changes in body weight and/or activity level.

### Cleaning

Use a damp cloth and mild soap to clean outside surfaces, do not use aggressive cleansers.

***The remaining instructions in this section are for practitioner use only.***

Maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

Carry out the following routine maintenance at least annually:

- Ensure all screws are secure. If not remove and clean screws and tighten screw to the correct torque settings.
- Check for visual defects that may affect proper function.
- Check for connection tightness. If loose, remove, clean and re-fit.

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check is recommended and signs of wear that may affect function should be reported to their service provider.

Advise the user to inform the practitioner/service provider of any changes in body weight and/or activity level.

If this device is used for extreme activity, the maintenance level and interval should be reviewed and if required advice and technical support sought to plan a new maintenance schedule dependant upon the frequency and nature of the activity. This should be determined by a local risk assessment carried out by a suitably qualified individual.

---

## 6 Limitations on Use

### Intended Life

A local risk assessment should be carried out based upon activity and usage.

### Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

### Environment

Avoid exposing the device to corrosive elements such as water, acids and other liquids. Also avoid abrasive environment such as those containing sand for example as these may promote premature wear.

Exclusively for use between -15 °C and 50 °C  
(5 °F to 122 °F).


Suitable for outdoor use

## 7 Bench Alignment

*The instructions in this section are for practitioner use only.*

Align with alignment slots medially. For use as part of a prosthetic device align as close as possible to the build line to maximise alignment possibilities and reduce excess loading to the device.

## 8 Fitting Advice



*The instructions in this section are for practitioner use only.*

Symptom	Remedy
A recurring noise occurs between the device and devices connected to it.	Ensure locking screw and single bolt connection are securely tightened to the torque ratings as specified in this document. Apply Loctite.
The adapter moves out of position.	User must not use the device until adjusted, repaired or replaced.

### Fitting Length

If the attaching tube needs to be cut to length, follow the instructions for use supplied with the tube, paying attention to the safety instructions associated with cutting the particular material of the tube.

Ensure the tube is cut square and the end has sharp edges removed.

---

## 9 Technical Data

Operating and  
Storage Temperature Range:

-15°C to 50°C  
(5°F to 122°F)

Component Weight (*size 26*):

130g

Activity Level:

1-3

Maximum User Weight:

100 kg (220 lb)

Range of Adjustment:


±360° rotation  
±5 mm horizontal shift

Build Height:

25 mm

---

### Build Height

## 10 Ordering Information

Spare Parts	
BK Top housing 30 mm	339008

### Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

### CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: [www.blatchford.co.uk](http://www.blatchford.co.uk)

Medical Device

Single Patient – multiple use

### Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

### Warranty

This device is warranted for 24 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

The above may vary by market; consult your local representative for details.

See the Blatchford website for the current full warranty statement.

## **Reporting of Serious Incidents**

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

## **Environmental Aspects**

This product is constructed from recyclable materials. Where possible the components should be recycled in accordance with local waste handling regulations.

## **Retaining the Packaging Label**

The practitioner is advised to keep the packaging label as a record of the device supplied.

## **Trademark Acknowledgements**

Registered trademarks of Blatchford Products Limited.

## **Manufacturer's Registered Address**

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**Blatchford Europe GmbH**
**Endolite India Ltd.**