Elite^{VT}

Instructions for Use

EVT22L11S–EVT30R94S
EVT22L11SD–EVT30R94SD

| | | |
|----|-----------------------|-----|
| EN | Instructions for Use | 2 |
| DA | Brugsanvisning | 17 |
| NO | Bruksanvisning | 32 |
| FI | Käyttöohjeet | 47 |
| SV | Bruksanvisning | 62 |
| RU | Инструкция протезиста | 77 |
| ZH | 使用说明 | 92 |
| AR | إرشادات الاستخدام | 107 |

Contents

| | |
|--|----|
| Contents | 2 |
| 1 Description and Intended Purpose | 3 |
| 2 Safety Information | 5 |
| 3 Construction | 6 |
| 4 Function..... | 7 |
| 5 Maintenance | 7 |
| 6 Limitations on Use | 8 |
| 7 Bench Alignment..... | 9 |
| 7.1 Static Alignment..... | 9 |
| 7.2 Dynamic Alignment | 9 |
| 7.3 Transfemoral Alignment..... | 9 |
| 8 Fitting Advice..... | 10 |
| 9 Foot Shell Removal & Foot Spring Replacement | 11 |
| 10 Axial Shock Spring Replacement..... | 13 |
| 11 Technical Data | 14 |
| 12 Ordering Information | 15 |

1 Description and Intended Purpose

These Instructions for Use are intended for use by the practitioner and user unless otherwise stated.

The term *device* is used throughout these instructions for use to refer to EliteVT.

Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.

Application

This device is to be used exclusively as part of a lower limb prosthesis.

Intended for a single user.

A high-energy-return foot. Independent heel and toe spring provides axial deflection. The split toe provides good ground compliance. This device also includes a spring element capable of absorbing axial and rotational shock intended to reduce shear stresses at the residuum/socket interface.

Features

- Precision titanium spring for vertical compression and axial rotation
- Linear response to axial loading
- Efficient energy return
- Independent e-carbon foot springs
- Heel wedge for adjusting spring stiffness

Activity Level

This device is recommended for users that have the potential to achieve Activity Level 3 or 4. Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.


Contraindications

This device may not be suitable for Activity Level 1 individuals or for competitive sports events, as these types of users will be better served by a specially designed prosthesis optimized for their needs.

Clinical Benefits

- Allows variable running speeds
- Increased walking speed, compared to rigid pylons
- High levels of mobility for active users, compared to non-energy-storage-and-return feet, hydraulic ankles, flex-foot-style feet and flex-walk-style feet.
- Improved prosthetic push-off work compared to solid-ankle-cushioned-heel feet
- Increased prosthetic positive work done, compared to lower activity feet
- High degree of user satisfaction, particularly with high activity users
- Reduced loading rate on the residual limb, compared to rigid pylons
- Reduced back pain when twisting, compared to rigid pylons

Spring Set Selection

- Low** Daily walking and occasional sports such as golf and hiking
- Moderate** Aggressive walking, frequent or daily sports such as jogging
- High** Daily activities such as distance running, climbing, lifting and carrying heavy objects for vocational purposes

Warning: For higher impact users, do not exceed the weight limit for individual springs.

Note... If in doubt choosing between two categories, choose the higher rate spring set.

Foot spring set recommendations shown are for transtibial users.

For transfemoral users we suggest selecting a foot spring set one category lower.

Once fitted with springs, cover appropriate lines on carrier with permanent black marker to leave spring set number showing.

2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the limb e.g. restricted or excessive movement, non-smooth motion or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

The device is not suitable for extreme sports, competitive running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk. Recreational cycling or running is acceptable.

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified practitioner.

The user should be advised to contact their practitioner if their condition changes.

Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.

To minimize the risk of slipping and tripping, appropriate footwear that fits securely onto the foot shell must be used at all times.

The device is not intended for use when immersed in water or as a shower prosthesis. If the limb comes into contact with water wipe dry immediately.

Be aware of finger trap hazard at all times

4 Function

This device comprises an e-carbon toe and independent heel spring. The carrier assembly contains a spring component which allows the pyramid tube to deflect and/or rotate through a limited range under an axial and/or torsional load. When the load is removed the spring returns to its original position. The effect of this movement is to attenuate any shock forces that would otherwise be transmitted to the residuum/socket interface both in axial and rotational directions.

Heel and toe springs are attached to the carrier assembly using titanium screws. The foot is wrapped in a UHM PE sock which is in turn surrounded by a PU foot shell.

5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider e.g. unusual noises, increased stiffness or restricted/excessive rotation, significant wear or excessive discoloration from long term exposure to UV.

Inform the practitioner/service provider of any changes in body weight and/or activity level.

Cleaning

Use a damp cloth and mild soap to clean outside surfaces, do not use aggressive cleansers.

The remaining instructions in this section are for practitioner use only.

This maintenance must be carried out only by competent personnel (practitioner or suitable trained technician).

The following routine maintenance is to be carried out at least annually:

- Remove foot shell and glide sock, check for damage or wear and replace if necessary.
- Check all screws for tightness, see *Construction* section; clean and reassemble as necessary.
- Check heel and toe springs for signs of delamination or wear and replace if necessary. Some minor surface damage may occur after a period of use, this does not affect the function or strength of the foot.

Instructions to remove the foot shell and replace the foot-springs are included in Foot Shell Removal & Foot Spring Replacement section.

If necessary, re-lubricate the spring's bearing surfaces. To do this:

- Disassemble the device as shown in the Axial Shock Spring Replacement instructions.
- Apply 'Sapphire Endure' (928015) or equivalent PTFE loaded grease to the bearing surface, then reassemble.

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check of the foot is recommended and signs of wear that may affect function should be reported to their service provider (e.g. significant wear or excessive discoloration from long term exposure to UV).

If this device is used for extreme activity, the maintenance level and interval should be reviewed and if required advice and technical support sought to plan a new maintenance schedule dependent upon the frequency and nature of the activity. This should be determined by a local risk assessment carried out by a suitably qualified individual.

6 Limitations on Use

Intended Life

A local risk assessment should be carried out based upon activity and usage.

Lifting Loads

Users weight and activity is governed by the stated limits.

Load carrying by the users should be based on a local risk assessment.

Environment

Avoid exposing this device to corrosive elements such as water, acids and other liquids. Also avoid abrasive environments such as those containing sand for example as these may promote premature wear.

Exclusively for use between -15 °C and 50 °C
(5 °F and 122 °F).


Suitable for outdoor use

7 Bench Alignment

The instructions in this section are for practitioner use only.

7.1 Static Alignment

Set-up Length

With flexion, adduction and abduction properly accommodated, set the limb length 10 mm longer than the sound side to allow for compression and deflection of the foot & axial shock springs during gait. This should be re-assessed before the dynamic trial commences and the length adjusted accordingly.

Build Line

This should fall 10 mm anterior to the centre of the pyramid (with heel height properly accommodated). The socket should be positioned accordingly.

7.2 Dynamic Alignment

Coronal Plane

Ensure that M-L thrust is minimal by adjusting relative positions of socket and foot.

Sagittal Plane

Check for smooth transition from heel strike to toe-off. Ensure also that when standing the heel and toe are evenly loaded and that both are touching the floor.

7.3 Transfemoral Alignment

Align transfemoral components according to fitting instructions supplied with the knee, keeping the build line relative to this device, as shown.

#Allow for users own footwear

*Approximate Ratio

8 Fitting Advice

The instructions in this section are for practitioner use only.

Springs are supplied as matched sets i.e. the heel, toe and axial springs are designed to work together to give smooth progression for most users.

Heel Wedge

A heel wedge is supplied with the foot. Fitting the wedge will have the effect of stiffening the heel spring. These can be taped in place for trial. For permanent fitting, wedges should be adhered in place by application of Loctite 424 (926104) between the lower contacting surface of the heel and the wedge.

Heel Stiffness

Progression throughout the stance phase should be smooth; heel function is key to this process:

- Too soft a heel or load line excessively posterior will result in sinking at heel strike and difficulty in getting over the toe.
- Too hard a heel or load line excessively anterior will result in a rapid progression through mid-stance or jarring at heel strike.

Axial Stiffness

The axial spring controls both compressive and torsional resilience. Typical vertical movement should be between 3-6 mm for normal walking. Alternative stiffness springs are available if required.

| | Symptoms | Remedy |
|----------------------------|--|---|
| Heel too soft | <ul style="list-style-type: none">• Sinking at heel strike• Difficulty climbing over the toe (toe feels too hard) | <ol style="list-style-type: none">1. Add heel wedge2. Move socket anteriorly in relation to the foot (excess movement may result in drop-off)3. If 1 and 2 fail, fit a stiffer spring set |
| Heel too hard | <ul style="list-style-type: none">• Rapid transition from heel strike through stance phase• Difficulty in controlling heel action, foot jars into mid-stance• Foot feels too rigid | <ol style="list-style-type: none">1. Remove heel wedge (if fitted)2. Move socket posteriorly in relation to foot3. If 1 and 2 fail, fit softer spring set |
| Toe/ Blade too soft | <ul style="list-style-type: none">• Rapid progression through mid-stance.• 'Drop off' at higher Activity Levels | <ol style="list-style-type: none">1. Move socket posteriorly in relation to foot2. Plantarflex foot slightly - note: some re-alignment may be required |

Note... Please contact your supplier if it is not possible to achieve a smooth gait after following the advice above.

9 Foot Shell Removal & Foot Spring Replacement

The instructions in this section are for practitioner use only.

Use appropriate health and safety equipment at all times including extraction facilities.

Be aware of finger trap hazard at all times.

1

13A/F

Remove toe spring screws (13A/F spanner 940273).

2


Withdraw toe spring from shell.

3

Rotate the carrier/heel spring assembly toward the rear of the foot to dis-engage the spring from its location in the shell.

4

Undo and remove heel spring screws. Assemble new heel spring onto carrier, use Loctite 243 (926012) and torque to 15 Nm. Use special Allen wrench 940080, Torque spanner adaptor 940081 or 13A/F spanner 940273

5

Attach the toe spring to the carrier using Loctite 243 (926012) on the bolts.

1. For spring rates 1 to 7 use 4 A/F Allen wrench and torque to 15 Nm. Do not use external hex, this is reserved for loosening the bolt, if required.

2. For spring rates 8 and 9 use 13 A/F spanner and torque to 25 Nm.

6

If a foam cosmesis is to be fitted, roughen top surface of foot shell to provide ideal bonding surface.

7

Mark carrier with spring category as appropriate and fit glide sock as shown.


8

Insert assembly into foot shell, push toe spring, carrier and heel assembly forward in the foot shell as far as possible.

9

Toe Spring location
in Foot Shell

Ensure toe spring's correct fit into location slot in foot shell.

10


Slide a metal shoe horn between heel spring and back of foot shell opening. Lever the spring down into the foot shell while pushing the whole assembly downwards.

11

The shoe horn should be located approx. 30 mm below top surface of foot shell to ensure full engagement of heel spring.

12

Ensure heel spring's correct fit into location slot in foot shell.

13





Bond foam cosmesis to top surface of the foot shell as shown using Thixofix adhesive (926204) or equivalent and shape to suit.

10 Axial Shock Spring Replacement

The instructions in this section are for practitioner use only.

Remove from foot shell as shown previously


Re-lubricate bearing surface with 'Sapphire Endure' (928015) or equivalent PTFE loaded grease if necessary.

Re-assembly is the reverse of above, note the position and location of anti-rotation pins in the axial spring and the orientation of the slot in the pylon (posterior). Apply Loctite and apply torque to screws as shown in the *Construction* section.

11 Technical Data

| | |
|--|-------------------------------------|
| Operating and Storage Temperature Range: | -15 °C to 50 °C (5 °F to 122 °F) |
| Component Weight (size 26): | 1000 g (2 lb 2 oz) |
| Activity Level: | 3-4 |
| Maximum User Weight: | 166 kg (365 lb) |
| Range of Adjustment | ±7° angular |
| Maximum Rotation | ±15° |
| Typical Axial Movement | 3-6 mm (max. 10 mm) |
| Proximal Attachment: | Male pyramid (Blatchford) |
| Build Height: | 170 mm |
| Heel Height | 10 mm |

Fitting Length

Note...
Refer to *Bench Alignment* section.

- * Sizes
- 22-26 = 65 mm
- 27-28 = 70 mm
- 29-30 = 75 mm

12 Ordering Information

Order Example

| | | | | | | |
|------------|-----------|---------------|-----------------|---------------|-----------------|---------------|
| EVT | 25 | L | N | 7 | 3 | S |
| | Size | Side (L/R) | Width* (N/W) | Spring Set | Shock Spring | Sandal Toe |

Available from
size 22 to size 30:
EVT22L11S to EVT30R94S
EVT22L11SD to EVT30R94SD
(add 'D' for a dark tone foot shell)

e.g. EVT25LN73S

*Sizes 25-28 only. For all other sizes, omit the Width field.

Spring Kit

| Spring Kit | Foot Size | | | |
|--------------|-----------|---------|---------|---------------|
| | 22-24 | 25-26 | 27-28 | 29-30 |
| Set 1 | 539701S | 539710S | 539719S | Special order |
| Set 2 | 539702S | 539711S | 539720S | Special order |
| Set 3 | 539703S | 539712S | 539721S | 539730S |
| Set 4 | 539704S | 539713S | 539722S | 539731S |
| Set 5 | 539705S | 539714S | 539723S | 539732S |
| Set 6 | 539706S | 539715S | 539724S | 539733S |
| Set 7 | 539707S | 539716S | 539725S | 539734S |
| Set 8 | 539708S | 539717S | 539726S | 539735S |
| Set 9 | | 539718S | 539727S | 539736S |

Foot Shell (For dark add 'D')

| Size/Side | Narrow | Wide |
|-----------|----------|----------|
| 22L | 539038S | - |
| 22R | 539039S | - |
| 23L | 539040S | - |
| 23R | 539041S | - |
| 24L | 539042S | - |
| 24R | 539043S | - |
| 25L | 539044SN | 539044SW |
| 25R | 539045SN | 539045SW |
| 26L | 539046SN | 539046SW |
| 26R | 539047SN | 539047SW |
| 27L | 539048SN | 539048SW |
| 27R | 539049SN | 539049SW |
| 28L | 539050SN | 539050SW |
| 28R | 539051SN | 539051SW |
| 29L | - | 539052S |
| 29R | - | 539053S |
| 30L | - | 539054S |
| 30R | - | 539055S |

Axial Shock Springs

| | |
|---------------|--------|
| Rate 1 | 539058 |
| Rate 2 | 539059 |
| Rate 3 | 539060 |
| Rate 4 | 539061 |

Glide Sock

| | |
|--------------------|--------|
| Sizes 22-24 | 531011 |
| Sizes 25-30 | 532811 |

Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: www.blatchford.co.uk

Medical Device

Single Patient – multiple use

Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

Warranty

This device is warranted for 36 months - foot shell 12 months - glide sock 3 months. The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions. See Blatchford website for the current full warranty statement.

Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

Environmental Aspects

Where possible the components should be recycled in accordance with local waste handling regulations.

Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

Trademark Acknowledgements

Elite Foot and Blatchford are registered trademarks of Blatchford Products Limited.

Manufacturer's Registered Address

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Blatchford Europe GmbH
Endolite India Ltd.