# Child's 4-Bar Knee

## Instructions for Use

239133

EN	Instructions for Use	2
DA	Brugsanvisning	12
NO	Bruksanvisning	22
FI	Käyttöohjeet	32
SV	Bruksanvisning	42
RU	Инструкция протезиста	52
ZH	使用说明	62
AR	إرشادات الاستخدام	72

# Contents

Contents .....	2
1 Description and Intended Purpose .....	3
2 Safety Information .....	5
3 Construction .....	6
4 Function.....	6
5 Maintenance .....	7
6 Limitations on Use .....	7
7 Bench Alignment.....	8
7.1 Static Alignment.....	8
7.2 Dynamic Alignment .....	8
8 Fitting Advice.....	8
8.1 Changing the Spring.....	8
8.2 Side Play Adjustment in the Front and Rear Links .....	9
9 Technical Data .....	10
10 Ordering Information .....	10

---

# 1 Description and Intended Purpose

These Instructions for Use are intended for use by the practitioner and user unless otherwise stated.

The term *device* is used throughout these instructions for use to refer to a Child's 4-Bar Knee.

Please read and ensure you understand all instructions for use, in particular all safety information and maintenance instructions.

## Application

This device is to be used exclusively as part of a lower limb prosthesis.

This device is for a single user.

It is a low profile polycentric knee for children which uses a 4 bar mechanism to provide geometric stability during stance with an adjustable extension assist for control during swing phase, and suitable for transfemoral and knee disarticulation users.

The geometry of the device provides increased toe clearance during swing phase. The build height is short and upon sitting minimizes protrusion of the prosthetic knee beyond the contra-lateral limb. This combination of features reduces asymmetry for users with longer residual limbs.

## Features

- Adjustable spring extension assist
- Stable geometry on heel strike
- Lightweight and robust knee
- A range of proximal adapters

## Activity Level

This device is recommended for users that have Activity Level 1-4. (Weight limits apply, see *Technical Data*.)

Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

### Activity Level 1

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

### Activity Level 2

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

### Activity Level 3

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

### Activity Level 4

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

## **Clinical Benefits**

- Improved ground clearance reducing trips and falls risk, compared to single axis knees
- Geometry increases knee stability during stance phase, compared to single axis knees
- Meets the needs of pediatric patients
- Improved aesthetics for knee disarticulation users or long transfemoral users, compared to single axis knees (e.g. does not stick out when seated)

## **Contraindications**

This device is not suitable for users with:

- Foot sizes larger than 21cm.
- Residual muscular weakness, contractures that cannot be correctly accommodated, or proprioceptive dysfunction including poor balance
- Contralateral joint instabilities or pathology
- Complicated conditions involving multiple disabilities

Users must be given gait training before using this device.

---

## 2 Safety Information

 This warning symbol highlights important safety information.

-  Beware of finger trap hazard all times. The use of cosmetics is recommended to minimise the risk of finger trap.
-  Any changes in the performance or function of the limb e.g. instability, double-action, restricted movement, non-smooth motion or unusual noises should be immediately reported to your service provider.
-  Always use a hand rail when descending stairs and at any other time if available.
-  Any excessive changes in the heel height after finalization of alignment may adversely affect limb function.
-  The device is not intended for use when immersed in water or as a shower prosthesis. If the limb comes into contact with water wipe dry immediately.
-  Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.
-  Care should be taken when carrying heavy loads as this may adversely affect the stability of the device.
-  The user should be advised to contact their practitioner if their condition changes.
-  Avoid exposure to extreme heat and/or cold.
-  Assembly, maintenance and repair of the device must only be carried out by a suitably qualified practitioner.


---

## 3 Construction

### Principal Parts

- Knee Assembly Aluminum, stainless steel, polymer, rubber
- Spring Stainless steel
- Adjuster Nylon, stainless steel
- Tube Clamp Screw Steel
- Lock Washer Steel

### Component Identification

---

## 4 Function

The device is a durable and lightweight pediatric knee designed for children of activity levels 1 to 4. The knee has a polycentric design for stability during stance, comprising pivots connecting a chassis to a main body via 4 link arms. The geometry of the knee provides efficient progression into swing phase at the required point in the gait cycle, increased toe clearance during the swing phase when compared to a single axis knee and short protrusion when sitting.

There is a choice of three extension assist springs with adjustable pre-load to provide swing control. Main bearings in the linkage are both replaceable and adjustable to prolong service life.

---

## 5 Maintenance

Visually check the device regularly.

Report any changes in performance of this device to the practitioner/service provider e.g. unusual noises, increased stiffness or restricted/excessive rotation or significant wear.

Inform the practitioner/service provider of any changes in body weight and/or activity level.

### Cleaning

Use a damp cloth and mild soap to clean the outside surfaces. DO NOT use aggressive cleansers. Thoroughly dry before use.

***The remaining instructions in this section are for practitioner use only.***

Maintenance must be carried out by competent personnel.

The following routine maintenance is to be carried out at least annually:

- Check the alignment bolt and tube clamp screw for tightness, tighten to the correct torque setting if required
- Check for defects, including corrosion that could affect proper function

To adjust any side play in the 4-bar links see section 8.2 *Side Play Adjustment in the Front and Rear Links*.

Advise the user to report any changes in the performance of this device to the practitioner.

Changes may include: loss of alignment, unusual noises, stiffness or looseness/play

Ensure the user has read and understood all safety and user-level maintenance information.

Advise the user that a regular visual check of the foot is recommended and any changes in performance and signs of wear that may affect function should be reported to their service provider (e.g. significant wear).

Advise the user to inform the practitioner/service provider of any changes in body weight and/or activity level.

---

## 6 Limitations on Use

### Intended Life

A local risk assessment should be carried out based upon activity and usage.

### Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

### Environment

Avoid exposing the product to corrosive elements such as water, acids and other liquids. Also avoid abrasive environments such as those containing sand for example as these may promote premature wear.


Exclusively for use between -10 °C and 50 °C (14 °F to 122 °F).

Suitable for outdoor use

## 7 Bench Alignment

*The instructions in this section are for practitioner use only.*

### 7.1 Static Alignment

#### Build Line

Construct the socket to fully accommodate alignment of residuum so that the Build Line passes through the attachment point.

The Build Line should fall in the middle of the slot on top of the Child's Knee, 12 mm anterior of upper front pivot.

Before walking stand the wearer in suitable walking rails to confirm alignment before progressing to dynamic alignment.

### 7.2 Dynamic Alignment

#### Adjust the Extension Assist

The adjuster can be used to customize the pre-load on the spring to suit the user's requirement. Undo the adjuster to reduce the assistance up to a maximum of 8 mm.

If satisfactory swing control cannot be achieved with the spring supplied. Fit one of the alternatives supplied with the kit.


## 8 Fitting Advice

### 8.1 Changing the Spring

To change the spring:

Fully unscrew the adjuster and remove the spring. Insert the appropriate spring and screw the adjuster back in with 1.5 mm showing.

The adjuster can then be used to adjust the pre-load on the spring to suit the user's individual requirements.



Spring	Spring Force	Extension Assist
Silver	Standard (as fitted)	Medium assistance
Gold	Weaker spring force	Low assistance
Black	Stronger spring force	High assistance



## 8.2 Side Play Adjustment in the Front and Rear Links

After a period of use it may be necessary to adjust the bearings to reduce side play in the links.


**3** To reduce the play in the Links:  
Using a Peg Spanner (Part No. 940054) in the Adjuster Cap to prevent rotation, use a second Peg Spanner to carefully turn the Pivot Pin clockwise to eliminate unnecessary side play in the links.  
Check that the knee moves freely.  
Note... Adjust one axis at a time and confirm free movement of knee before progressing to adjust the next.  
Note... If play cannot be eliminated replace bearings and apply STP to lubricate bearing surfaces.

## 9 Technical Data

Operating and Storage Temperature Range:	-15 °C to 50 °C (5 °F to 122 °F)
Component Weight:	335 g (12 oz)
Activity Level:	1-4
Maximum User Weight:	60 kg (130 lb)
Distal Attachment:	Ø 25 mm Full 360° rotation
Proximal Attachment:	Blatchford Single Bolt with 7 mm A-P shift
Build Height:	105 mm
[See diagram below]	

### Fitting Length

\*Maximum protrusion inside shin.

## 10 Ordering Information

### Order Example

Item	Part Number
Child's 4-Bar Knee	239133

### Parts and Tools

Item	Part Number
Peg Spanner	940054
Bearing Kit	339006
Bearing and Links Kit	339010
Tube Clamp Screw	910215

Spring Kit		Part Number
Spring	Silver	925014
	Gold	925015
	Black	925016

## Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

## CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: [www.blatchford.co.uk](http://www.blatchford.co.uk)

Medical Device

Single Patient – multiple use

## Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a practitioner.

## Warranty

This device is warranted for 24 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licenses and exemptions.

See the Blatchford website for the current full warranty statement.

## Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

## Environmental Aspects

Where possible the components should be recycled in accordance with local waste handling regulations.

## Retaining the Packaging Label

The practitioner is advised to keep the packaging label as a record of the device supplied.

## Manufacturer's Registered Address

Blatchford Products Limited, Lister Road, Basingstoke RG22 4AH, UK.

## Trademark Acknowledgements

Blatchford is a registered trademark of Blatchford Products Limited.

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**Blatchford Europe GmbH**
**Endolite India Ltd.**