# Elite<sup>2</sup>

## Instructions for Use

EL22L1S–EL30R9S  
EL22L1SD–EL30R9SD

EN	Instructions for Use	2
SL	Navodila za uporabo	16
BG	Инструкции за употреба	30
HR	Upute za uporabu	44
SK	Návod na použitie	58
HU	Használati útmutató	72
EL	Οδηγίες χρήσης	86
LV	Lietošanas pamācība	100
LT	Naudojimo instrukcija	114
ET	Kasutusjuhend	128

# Contents

Contents .....	2
1 Description and Intended Purpose .....	3
2 Safety Information .....	5
3 Construction .....	6
4 Function.....	7
5 Maintenance .....	7
6 Limitations on Use .....	7
7 Bench Alignment.....	8
8 Fitting Advice.....	9
9 Dis-assembly Instructions .....	10
10 Assembly Instructions .....	11
11 Technical Data .....	13
12 Ordering Information .....	14

---

# 1 Description and Intended Purpose

These instructions for use are for the practitioner.

The term *device* is used throughout these instructions for use to refer to Elite2.

## **Application**

This device is to be used exclusively as part of a lower limb prosthesis.

Intended for a single user.

A high-energy-return foot. The independent heel and toe springs provide axial deflection. The split toe provides good ground compliance.

## **Activity Level**

This device is recommended for users that have the potential to achieve Activity Level 3 or 4. Of course there are exceptions and in our recommendation we want to allow for unique, individual circumstances and any such decision should be made with sound and thorough justification.

### **Activity Level 1**

Has the ability or potential to use a prosthesis for transfers or ambulation on level surfaces at fixed cadence. Typical of the limited and unlimited household ambulator.

### **Activity Level 2**

Has the ability or potential for ambulation with the ability to traverse low-level environmental barriers such as curbs, stairs, or uneven surfaces. Typical of the limited community ambulator.

### **Activity Level 3**

Has the ability or potential for ambulation with variable cadence. Typical of the community ambulator who has the ability to traverse most environmental barriers and may have vocational, therapeutic, or exercise activity that demands prosthetic utilization beyond simple locomotion.

### **Activity Level 4**

Has the ability or potential for prosthetic ambulation that exceeds basic ambulation skills, exhibiting high impact, stress, or energy levels. Typical of the prosthetic demands of the child, active adult, or athlete.

## Contraindications

This device might not be suitable for Activity Level 1 individuals or for competitive sports events, as these types of users will be better served by a specially designed prosthesis optimized for their needs.

Ensure that the user has understood all instructions for use, drawing particular attention to the section regarding maintenance.

## Clinical Benefits

- Vertical compliance for shock-absorption
- High levels of mobility for active users
- Improved prosthetic side push-off
- Increased walking speed

## Spring Set Selection

	User Weight [kg (lbs)]									
	44-52 (100-115)	53-59 (116-130)	60-68 (131-150)	69-77 (151-170)	78-88 (171-195)	89-100 (196-220)	101-116 (221-255)	117-130 (256-285)	131-147 (286-325)	148-166 (326-365)
Impact	Activity Level 3 : Foot Spring Set									
Low	1	1	2	3	4	5	6	7	8	9
	Activity Level 4 : Foot Spring Set									
Mod	1	2	3	4	5	6	7	8	9	
High	2	3	4	5	6	7	8	9		

Low Daily walking and occasional sports such as golf and hiking

Moderate Aggressive walking, frequent or daily sports such as jogging

High Daily activities such as distance running, climbing, lifting and carrying heavy objects for vocational purposes

### Important:

**For higher impact users, do not exceed the weight limit for individual springs.**

Note... If in doubt choosing between two categories, choose the higher rate spring set.

**Foot Spring set recommendations shown are for transtibial users.**

**For transfemoral users we suggest selecting a spring set one category lower, refer to Section 8 *Fitting Advice*, to ensure satisfactory function and range of movement**

---

## 2 Safety Information

This warning symbol highlights important safety information which must be followed carefully.

Any changes in the performance or function of the limb e.g. restricted movement, non-smooth motion or unusual noises should be immediately reported to your service provider.

Always use a hand rail when descending stairs and at any other time if available.

The device is not suitable for extreme sports, running or cycle racing, ice and snow sports, extreme slopes and steps. Any such activities undertaken are done so completely at the users' own risk. Recreational cycling is acceptable.

Assembly, maintenance and repair of the device must only be carried out by a suitably qualified clinician.

The user should contact their clinician if their condition changes.

Ensure only suitably retrofitted vehicles are used when driving. All persons are required to observe their respective driving laws when operating motor vehicles.


To minimise the risk of slipping and tripping, appropriate footwear that fits securely onto the footshell must be used at all times.

# 3 Construction

## Principal Parts:

- Carrier Assembly (aluminum/stainless steel/titanium)
- Heel & Toe Springs (e-Carbon)
- Spring Attachment Screws (titanium)
- Glide Sock (UHM PE)
- Foot Shell (PU)

---

## 4 Function

This device comprises an e-carbon toe and independent heel spring.

Heel and toe springs are attached to the carrier using titanium screws. The foot is wrapped in a UHM PE sock which is in turn surrounded by a PU foot shell.

---

## 5 Maintenance

Maintenance must be carried out by competent personnel.

We recommend the following annual maintenance:

- Remove the foot shell and glide sock, check for damage or wear and replace if necessary.
- Check all screws for tightness, clean and reassemble as necessary.
- Check heel and toe springs for signs of delamination or wear and replace if necessary. Some minor surface damage may occur after a period of use, this does not affect the function or strength of the foot.

The user should be advised:

- Any changes in performance of this device must be reported to the practitioner e.g. reduced energy return or unusual noises.
- The practitioner must also be informed of any changes in body weight and/or activity level.

If this device is used for extreme activity, the maintenance level and interval should be reviewed and if required advice and technical support sought to plan a new maintenance schedule dependent upon the frequency and nature of the activity. This should be determined by a local risk assessment carried out by a suitably qualified individual.

The user should be advised that a regular visual check of the foot is recommended, signs of wear that may affect function should be reported to their service provider (e.g. significant wear or excessive discoloration from long term exposure to UV).

### Cleaning

Use a damp cloth and mild soap to clean outside surfaces, do not use aggressive cleansers.

---

## 6 Limitations on Use

### Intended Life

A local risk assessment should be carried out based upon activity and usage.

### Lifting Loads

User weight and activity is governed by the stated limits.

Load carrying by the user should be based on a local risk assessment.

### Environment

This device is waterproof to a maximum depth of 1 meter.

Thoroughly rinse this device with fresh water after use in abrasive environments such as those that may contain sand or grit, for example, to prevent wear or damage to moving parts.

Thoroughly rinse with fresh water after use in salt or chlorinated water.

Foot units must be adequately finished to prevent water ingress into the foot shell where possible. If water enters the foot shell, it should be inverted and dried before further use.


Exclusively for use between -15 °C and 50 °C.

Suitable for submersion

# 7 Bench Alignment

## 7.1 Static Alignment

### Setup length

With flexion, adduction and abduction properly accommodated, set the limb length 5 mm longer than the sound side to allow for compression and deflection of the foot springs during gait. This should be re-assessed once the dynamic trial commences and the length adjusted accordingly.

### Build Line

This should fall 10 mm anterior to the centre of the pyramid (with heel height properly accommodated). The socket should be positioned accordingly.

### Dynamic Alignment

#### Coronal Plane

Ensure that M-L thrust is minimal by adjusting relative positions of socket and foot.

#### Sagittal Plane

Check for smooth transition from heel strike to toe-off. Ensure also that when standing the heel and toe are evenly loaded and that both are touching the floor.

#### Transfemoral Alignment

Align trans-femoral components according to fitting instructions supplied with the knee, keeping the build line relative to the foot as shown.

# Allow for users own footwear

\* Approximate ratio

## 8 Fitting Advice

Spring sets are supplied as matched pairs i.e. the heel and toe spring are designed to work together to give smooth progression for most users.

### Heel Wedge

A heel wedge is supplied with the foot. Fitting the wedge will have the effect of stiffening the heel spring. These can be taped in place for trial. For permanent fitting, wedges should be adhered in place by application of Loctite 424 (926104) between the lower contacting surface of the heel and the wedge.

### Heel Stiffness

Progression throughout the stance phase should be smooth; heel function is key to this process:

- Too soft a heel or load line excessively posterior will result in sinking at heel strike and difficulty in getting over the toe.
- Too hard a heel or load line excessively anterior will result in a rapid progression through mid-stance or jarring at heel strike.

	Symptoms	Remedy
<b>Heel too soft</b>	<ul style="list-style-type: none"><li>• Sinking at heel strike</li><li>• Difficulty climbing over the toe (toe feels too hard)</li></ul>	<ol style="list-style-type: none"><li>1. Add heel wedge</li><li>2. Move socket anteriorly in relation to the foot (excess movement may result in drop off)</li><li>3. If 1 and 2 fail, fit a stiffer spring set</li></ol>
<b>Heel too hard</b>	<ul style="list-style-type: none"><li>• Rapid transition from heel strike through stance phase</li><li>• Difficulty in controlling heel action, foot jars into mid-stance</li><li>• Foot feels too rigid</li></ul>	<ol style="list-style-type: none"><li>1. Remove heel wedge (if fitted)</li><li>2. Move socket posteriorly in relation to foot</li><li>3. If 1 and 2 fail, fit softer spring set</li></ol>
<b>Toe/ Blade too soft</b>	<ul style="list-style-type: none"><li>• Rapid progression through mid-stance.</li><li>• 'Drop off' at higher Activity Levels</li></ul>	<ol style="list-style-type: none"><li>1. Move socket posteriorly in relation to foot</li><li>2. Plantarflex foot slightly - note: some re-alignment may be required</li></ol>

Please contact your supplier if it is not possible to achieve a smooth gait after following the advice above.

## 9 Dis-assembly Instructions

Use appropriate health and safety equipment at all times including extraction facilities.

Be aware of finger trap hazard at all times.

1

Carefully remove, with a knife, any foam cosmetics that may be adhered to the foot shell.

2


Remove securing screws from the carrier and pull the toe spring towards the rear of the foot.

3

Remove the toe spring to leave the carrier/heel spring assembly alone inside the foot shell.

4

Rotate the carrier/heel spring assembly toward the rear of the foot to dis-engage the spring from its location in the shell.

# 10 Assembly Instructions


**1**

Carrier  
Heel spring

Assemble heel spring onto carrier as shown.

**2**

LOCTITE 243

4  
15 Nm  
or  
13A/F

Use special Allen wrench 940080, Torque spanner adaptor 940081 or 13A/F spanner 940273. Use Loctite 243 (926012) and torque to 15 Nm.

**3**

LOCTITE 243

Attach the toe spring to the carrier using Loctite 243 (926012) on the bolts.

1. For spring rates 1 to 7 use 4 A/F Allen wrench and torque to 15 Nm. Do not use external hex, this is reserved for loosening the bolt, if required.

2. For spring rates 8 and 9 use 13 A/F spanner and torque to 25 Nm.

**4**

Ensure correct spring set number is marked on left side of carrier. Use black marker to cover unwanted segments.

**5**

Glide Sock

Fit glide sock onto the toe spring, heel spring and carrier assembly as shown.

**6**

If a foam cosmesis is to be fitted, roughen top surface of foot shell to provide ideal bonding surface.

**7**

Toe spring location in foot shell

Slide whole assembly into foot shell. Ensure correct fit into toe spring slot.

**8**

Slide a metal shoe horn between heel spring and back of foot shell opening and lever the spring down into the foot shell.

**9**

30 mm below heel

Press heel spring into location in shell as shown. The shoe horn should be located approx. 30 mm below top surface of foot shell

**10**

heel spring location slot

Ensure heel spring's correct fit into location slot in foot shell.


**11**

Bond foam cosmesis to top surface of the foot shell as shown using Thixofix adhesive (926204) or equivalent and shape to suit.

# 11 Technical Data

Operating and Storage Temperature Range:	-15 °C to 50 °C (5 °F to 122 °F)
Component Weight (size 26):	680 g (1 lb 7 oz)
Activity Level:	3–4
Maximum User Weight:	166 kg (365 lb)
Proximal Alignment Attachment:	Male Pyramid (Blatchford)
Range of Adjustment:	±7° Angular
Build Height: (See diagram below)	130 mm
Heel Height:	10 mm

## Fitting Length

\* Sizes  
22-26 = 65mm  
27-28 = 70mm  
29-30 = 75mm

# 12 Ordering Information

## Order Example

<b>EL</b>	<b>25</b>	<b>L</b>	<b>N</b>	<b>3</b>	<b>S</b>
	Size	Side (L/R)	Width* (N/W)	Spring Set Category	Sandal Toe

Available from size 22 to size 30:

EL22L1S to EL30R9S

EL22L1SD to EL30R9SD

(add 'D' for a dark tone foot shell)

e.g. EL25LN3S

\*Sizes 25–28 only. For all other sizes, omit the Width field.

## Spring Kit

Spring Kit	Foot Size			
	22-24	25-26	27-28	29-30
Set 1	539701S	539710S	539719S	Special order
Set 2	539702S	539711S	539720S	Special order
Set 3	539703S	539712S	539721S	539730S
Set 4	539704S	539713S	539722S	539731S
Set 5	539705S	539714S	539723S	539732S
Set 6	539706S	539715S	539724S	539733S
Set 7	539707S	539716S	539725S	539734S
Set 8	539708S	539717S	539726S	539735S
Set 9		539718S	539727S	539736S

## Foot Shell

(For dark add 'D')

Size/Side	Narrow	Wide
22L	539038S	-
22R	539039S	-
23L	539040S	-
23R	539041S	-
24L	539042S	-
24R	539043S	-
25L	539044SN	539044SW
25R	539045SN	539045SW
26L	539046SN	539046SW
26R	539047SN	539047SW
27L	539048SN	539048SW
27R	539049SN	539049SW
28L	539050SN	539050SW
28R	539051SN	539051SW
29L	-	539052S
29R	-	539053S
30L	-	539054S
30R	-	539055S

## Glide Sock

Sizes 22-24	531011
Sizes 25-30	532811

## Liability

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use supplied with the device. The manufacturer is not liable for any adverse outcome caused by any component combinations that were not authorized by them.

## CE Conformity

This product meets the requirements of the European Regulation EU 2017/745 for medical devices. This product has been classified as a class I device according to the classification rules outlined in Annex VIII of the regulation. The EU declaration of conformity certificate is available at the following internet address: [www.blatchford.co.uk](http://www.blatchford.co.uk)

Medical Device

Single Patient – multiple use

## Compatibility

Combination with Blatchford branded products is approved based on testing in accordance with relevant standards and the MDR including structural test, dimensional compatibility and monitored field performance.

Combination with alternative CE marked products must be carried out in view of a documented local risk assessment carried out by a Practitioner.

## Warranty

This device is warranted for 36 months - foot shell 12 months - glide sock 3 months.

The user should be aware that changes or modifications not expressly approved could void the warranty, operating licences and exemptions.

See Blatchford website for the current full warranty statement.

## Reporting of Serious Incidents

In the unlikely event of a serious incident occurring in relation to this device it should be reported to the manufacturer and your national competent authority.

## Environmental Aspects

Where possible the components should be recycled in accordance with local waste handling regulations.

## Retaining the Packaging Label

You are advised to keep the packaging label as a record of the device supplied.

## Manufacturer's Registered Address

Blatchford Products Limited, Lister Road, Basingstoke RG22 4AH.

## Trademark Acknowledgements

Elite Foot and Blatchford are registered trademarks of Blatchford Products Limited.

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**Blatchford Europe GmbH**
**Endolite India Ltd.**