

BACK TO SCHOOL!

Paediatric Footwear From Stock

Catapult Boot II



Sizes 27-45
Last Medium and Wide
Fastening Lace
Height UK
Stiffeners Standard
Colour Nubuck Black

£99.99

Part No. **Medium** 78021BK*M
Wide 78022BK*W * = Insert size required


Slickers Boot II



Sizes 27-45
Last Medium and Wide
Fastening Lace
Height UK
Stiffeners Standard
Colour Natural Leather Black

£99.99

Part No. **Medium** 78001BK*M
Wide 78002BK*W * = Insert size required



Sprinter Boot



Sizes 27-45
Last Medium and Wide
Fastening Lace
Height UK
Stiffeners Standard
Colour Nubuck Black

£109.99

Part No. **Medium** 78501BK*M
Wide 78502BK*W * = Insert size required

Boot Size	Medium					Wide				
	Length	Joint Width	Instep Circ.	Ankle Circ.	Toe Depth	Joint Width	Instep Circ.	Ankle Circ.	Toe Depth	
27	179	65	174	166	21	68	182	172	22	
28	186	67	178	171	22	70	186	177	22	
29	193	68	182	176	22	71	190	182	23	
30	199	69	186	180	23	72	194	188	23	
31	206	70	190	185	23	73	198	193	24	
32	213	72	194	190	24	75	202	198	25	
33	219	74	198	195	25	77	206	203	25	
34	226	75	202	200	25	78	210	208	26	
35	232	77	206	205	26	80	214	218	26	
36	239	79	210	210	26	82	218	224	27	
37	246	80	214	215	27	83	222	229	28	
38	252	82	218	220	28	85	226	235	28	
39	259	83	222	225	28	86	230	240	29	
40	266	85	226	230	29	88	234	246	29	
41	272	87	230	238	29	90	238	253	30	
42	279	88	234	246	30	91	242	259	31	
43	286	90	238	254	31	93	246	266	31	
44	292	92	242	252	31	95	250	272	32	
45	299	93	246	270	32	96	254	278	32	

Measurements for all sections are shown in mm.

All products will be shipped within 48 hours dependent on stock availability.
 All footwear is adaptable with a range of sole additions.

Blatchford

+44 (0) 114 2637900 | orthoticscs@blatchford.co.uk | blatchfordmobility.com | @blatchfordgrp